

## Device-Specific Handling Error Checklist (DSHEC).

Type	Spiromax®	Easyhaler®	Turbuhaler®
<b>Preparation</b>	Does not hold the device in up-right position with mouthpiece down	Does not remove the cap	Does not remove the cap
	Not fully opening the cap (a click is not heard)	Does not hold the device in up-right position with mouthpiece down	Holds the mouthpiece when priming (a click not heard)
		Does not shake the device up and down 3-5 times	Does not do proper rotation when priming the device (a click not heard)
		Priming the device at the same time as shaking it	Does not twist forward and back when priming
		A click is not heard when priming	Does not hold the device in up-right position after priming ( $\pm 45^\circ$ )
		Does not hold the device in up-right position when priming ( $\pm 45^\circ$ )	
		Holds the plunger down when inhaling	
	Does not empty the device when loading more than one dose		
	Does not empty the device when loading more than one dose		
<b>Inhalation</b>	Exhales into the device	Exhales into the device	Exhales into the device
	Fails to put in mouth and seal the lips around the mouthpiece	Fails to put in mouth and seal the lips around the mouthpiece	Fails to put in mouth and seal the lips around the mouthpiece
	Inhalation is not as fast/deep as possible from the beginning	Inhalation is not as fast/deep as possible from the beginning	Inhalation is not as fast/deep as possible from the beginning
	Inhales through the nose	Inhales through the nose	Inhales through the nose
	Covers the air inlets during the inhalation		
<b>After inhalation</b>	Does not place back the cap	Does not place back the cap	Does not place back the cap